


For more information:

Request a demonstration or seminar for
your lab

Carolyn@douglas.co.uk



Douglas Instruments
Success in protein crystallization

Dear Crystallographer

Use less protein... explore more crystallization space...


Save protein and increase drop accuracy and reproducibility by using an [Oryx robot](#) for screening and optimization. Because [almost no protein is wasted](#) by the robot more experiments can be dispensed and more crystallization space can be covered.

Use less seed stock

Seed stock is very valuable and for [rMMS screening experiments](#) it's best to use a high concentration of seed stock, so it's important not to waste seed stock.

Great for [scaling up](#) too

Drops are accurately dispensed using the Oryx's microtip, which can dispense volumes from 5nL** to 8 µL. Up to 5 drops can be dispensed to hanging drop cover slides allowing "mini-optimizations" in each well!



Crystal production and harvesting

For production of a large number of similar sized crystals, adding seed stock to the drop is an ideal way to improve reproducibility. There's an experiment script for [microseeding optimization](#) that finds the correct dilution of seed stock for up to 8 hit conditions on a single plate.

**such small volumes are only recommended for seeding experiments

Oryx protein usage

Method	Min drop volume	Max drop volume	Total protein required: 24 well optimization 0.2+0.2 µL	Total seed stock required: 24 wells, 20 nL
Hanging drop	100+100 nL	up to 8+8 µL	5 µL	0.5 µL
Sitting drop	100+100 nL	up to 8+8 µL	5 µL	0.5 µL
Microbatch-under-oil	100+100 nL	up to 8+8 µL	5 µL	0.5 µL

[Request more information](#)


Products available from Douglas Instruments



Oryx Robots



Crystallization Plates



Oryx Microtips

Conferences:

Douglas Instruments is attending the following meetings:

BCA spring meeting 2022, Leeds, UK
11 April - 14 April 2022

ACA 2022, Portland, USA
30 July - 2 August 2022

Recently published research using Oryx protein crystallization robots:

Microbial enzymes induce colitis by reactivating triclosan in the mouse gastrointestinal tract

Zhang, Jianan, et al., 2022.

Nature Communications, 13(1), pp.1-14.

Elucidating the 3D Structure of a Surface Membrane Antigen from *Trypanosoma cruzi* as a Serodiagnostic Biomarker of Chagas Disease

Di Pisa, Flavio, et al., 2022.

Vaccines, 10(1), p.71.

