


Lipidic Cubic Phase Mixer



- Manually operated device for mixing Lipidic Cubic Phase (LCP)
- Easy to use – create LCP without leaks or cracked syringes.
- Manual operation allows user to feel mixing process. Control speed of mixing and choose where to mix.
- Mix volumes of LCP from 10 μ l or less to 100 μ l.
- Minimal dead volume – No material is wasted in a mixer union.
- The LCP mixer kit is available separately or included with Oryx LCP option.
- Included in kit:

Mixer device; 2 Hamilton syringes; LCP dispensing needle (6.5mm polished); spare syringe fitting.